



Sonomed Escalon™

VuPad™

An innovation in ultrasound
you can see and touch



OPERATOR'S MANUAL

CONTENTS

OVERVIEW	4
<hr/>	
Description	
Indications for Use	
Contraindications	
GETTING STARTED	6
<hr/>	
System and Components	
Unpacking	
Consideration When Joining Network Domain	
Opening Display and Set-Up	
A-SCAN	10
<hr/>	
Patient Database	
Exam Mode	
Verify Calibration	
Exam Mode	
Calculation Mode	
Reports	
DIAGNOSTIC A-SCAN	14
<hr/>	
PACHYMETRY	17
<hr/>	
Select Exam Mode	
Perform Exam	
Report	
B-SCAN AND UBM	20
<hr/>	
Configure System	
Patient Database	
Exam Mode	
MAINTENANCE	23
<hr/>	
Cleaning and Disinfection	
SPECIFICATIONS / PERFORMANCE CHARACTERISTICS	26
<hr/>	
General Specifications	
A-Scan Specifications	
B-Scan Specifications	
Pachymeter Specifications	
UBM Specifications	
Data and Image Processing Specifications	
ALARA and Emissions	
Immunity Test Levels	
WARNINGS AND CAUTIONS	38
<hr/>	
CYBER SECURITY RECOMMENDATIONS	40
<hr/>	
WARRANTY	41
<hr/>	
SYMBOLS	42
<hr/>	

CONTACT INFORMATION



Sonomed, Inc.
1979 Marcus Avenue, C105
Lake Success, NY 11042, U.S.A.
Tel: 800-227-1285
Tel: 516-354-0900
Fax: 516-354-5902
info@sonomedinc.com
www.sonomedescalon.com



Authorized Representative for Medical Device Directive:
Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands



United States Federal law and European regulations restrict the use of this medical device to, or on the order of, a physician. The VuPad ophthalmic ultrasound device is for use by trained professionals in optometry and ophthalmology in a clinical setting. This group includes ophthalmologists, optometrists and ophthalmic technicians. The device is not intended for use by lay persons in any capacity.

A-scan biometry, also referred to as A-scan, utilizes an ultrasound device for diagnostic testing and is used as co-management tool in examinations of the structural integrity and pathology of the eye. The A-scan can determine the length of the eye and can be useful in diagnosing common sight disorders and is beneficial in cataract surgeries, as it can determine the power of the intraocular lens (IOL) needed for the artificial implant. In addition to axial length, ultrasound biometry can measure anterior chamber depth and lens thickness. Pachymetry is for measuring and mapping corneal thickness of the eye. The B-Scan mode produces a live, two-dimensional image to facilitate the identification and measurement of ocular pathologies in the posterior chamber of the eye, particularly when view of the chamber is obscured, such as is the case with cataracts. A UBM (ultrasound bio-microscope) is a B-scan which utilizes higher frequency transducers than used in a normal B-scan. Using higher frequencies provides much higher image resolution but limits the depth penetration of the signal and is therefore used primarily to visualize the internal components of the anterior segment of the eye only (iris, lens, etc.).

There are no restrictions to the patient population with regards to gender, age or physical limitations. The device is not for use on patients with questionable ocular integrity.

Before examining a patient, the user should become acquainted with the operating procedures, warnings and precautions set forth in the operator's manual. The user should consult additional resources as necessary for further information regarding the proper application of ultrasound technology. This instrument should be used in strict accordance with the instructions outlined in this operator's manual. The safety of the operator and the performance of the instrument cannot be guaranteed if used in a manner not specified by Sonomed Escalon.

Do not use the device together with HF surgical equipment. HF surgical equipment may be damaged, which may result in fire.

If difficulty is experienced when operating the unit after carefully reviewing this operator's manual, contact your local Sonomed Escalon distributor for assistance.

There are no user-serviceable parts within the VuPad system.



To receive a translated copy of this manual, contact your in-country distributor, or call Sonomed directly at 516-354-0900 or 800-227-1285. For technical service and support please contact Sonomed Escalon or your local distributor.

OVERVIEW

DESCRIPTION

The VuPad™ is an ophthalmic ultrasound system capable of utilizing a range of transducer frequencies for optimized imaging and measurements of various ophthalmic structures. The device is a multi-purpose system intended for use in ophthalmic applications designed to capture images of the interior of the eye and make accurate measurements of the structures.

The VuPad™ is a stand-alone system that runs on a Windows 10 platform and may be networked (by the user) for interface with electronic medical records systems, printing, and other purposes. The system consists of the VuPad™ console, ultrasound probe(s) and transducer(s), and foot pedal.

The device is used by coupling the probe/transducer to the eye either through direct contact or immersion methods. Available modes are biometric A-scan, pachymeter, diagnostic A-scan, B-scan, UBM (ultrasound bio-microscope).

Biometric A-Scan

A-Scan is ultrasound technology that is used to produce length measurements along the visual axis of the eye for calculation of intraocular lens power. This method enables precise measurement of the anterior chamber depth (ACD), lens (L), and vitreous to produce the axial length of the eye. When a cataract is removed, the lens is replaced with an artificial lens implant. By measuring both the axial length of the eye (A-scan) and the power of the cornea (keratometry), a user selected formula can be used to calculate the power of the intraocular lens needed.

Pachymeter

Pachymetry is ultrasound technology that is used to produce measurements of the thickness of the eye's cornea. It is used to perform corneal pachymetry prior to refractive surgery, for keratoconus screening, limbal relaxing incisions surgery and is useful in screening for patients suspected of developing glaucoma among other uses.

Diagnostic A-Scan

Diagnostic A-Scan is ultrasound technology that is used for diagnostic evaluation of detected eye pathologies found with B-Scan screening. This method enables imaging of structural amplitudes for analytical determination of the patient's eye disorder.

B-Scan

B-Scan is ultrasound technology that produces a cross-section, two-dimensional grayscale images for diagnosing pathologies of the posterior segment of the eye. This method enables imaging when the light-conducting media of the eye are opaque. Common conditions such as cataract, vitreous degeneration, retinal detachment, ocular trauma, choroidal melanoma, and retinoblastoma can be accurately evaluated with this modality.

UBM/ Ultrasound Bio-Microscopy (UBM-Mode)


UBM is ultrasound technology that produces high resolution cross-section, two-dimensional grayscale images of the anterior segment of the eye. This method enables imaging of structural details such as Bowman's membrane, stroma, cornea, anterior chamber, lens, iris, ciliary body, and scleral spur. Other parameters such as iris area and volume, angle opening distance, angle recess area, sclera thickness, and trabecular meshwork-ciliary process can be evaluated and measured. Common conditions such as glaucoma, iris cyst, neoplasms, trauma and foreign bodies can be accurately identified.

INDICATION FOR USE

The VuPad™ is intended to be used to visualize and measure the eye and orbit using A-Scan, B-Scan and pachymeter ultrasound.


CONTRAINDICATIONS

The VuPad™ is not intended for fetal use.



GETTING STARTED

SYSTEM AND COMPONENTS



Power On/Off

Ethernet

USB

Power Cable

A-Probe 10 MHz



B-Probe
(12 or 20 MHz Probes)



Water Path Probe
(35 or 50 MHz Transducers)



Pachymetry Probe 20 MHz



Diagnostic A-Probe 8 MHz



Foot Pedal



UNPACKING

- ① Unpack contents from packaging.
- ② Examine the unit and components to ensure the contents are intact. Visually examine the probes for any signs of cracks, scratches or damage. Do not use probes if damage is apparent.
- ③ Connect power cord between console and wall outlet.
NOTE: Console may be propped on its stand or attached to a VESA mount.
CAUTION: Position such that console is well ventilated with easy access to disconnect power cords, as may be necessary.
- ④ Connect the probe cables to unit and place probes into probe holders.
NOTE: The same cable is used for all B-scan and UBM probes.
- ⑤ Power on console. System will boot up into a Windows 10 home screen and the VuPad icon can be selected to launch the program.
NOTE: Console may be powered off by pushing the power button.

CONSIDERATION WHEN JOINING NETWORK DOMAIN

The Vupad uses Windows 10 and may be joined to (or removed from) a network domain. However, when doing so, it is required that the local user account is set up as a local administrator on the ultrasound system.

One potential issue to keep in mind is that when an ultrasound system is joined to a domain, the domain rules are typically pushed down, which can prevent a local user account from having full permissions (this would result in the ultrasound system not functioning properly). Please ensure that the local user account is set up as an administrator on the ultrasound system via your domain rules (i.e. the local user must remain as the local administrator).

For questions, please contact technical service at 516-354-0900 or 800-227-1285 or email ultrasound-support@escalonmed.com.




Fig. 2

Lens Selection Screen:

Manufacturer: Scroll and Select the Lens Manufacturer's Name.

Lenses: Scroll and Select the desired Lens Model.

NOTE: The "Lens Database" Tab (see Fig.1) includes over 1500 lenses and associated data which the user can update, by adding, deleting or editing lenses.

The "LENSES" Tab is the Examiner Lens Selection tab from which the Examiner must choose the desired lenses from for performing IOL calculations.

A-SCAN

PATIENT DATABASE

- ① **PATIENT:** Select the “NEW” tab located at the top left side of the database display screen.
NOTE: For previously entered patient, scroll through patient list, by sliding the scroll indicator just to the right of patient list until the desired patient name is displayed.
- ② **ENTER NEW PATIENT** Information into the appropriate fields (see Fig. 3)
[Note: K-readings can be entered here or later if performing a-scan biometry.]..... Press “OK”




Fig. 3

EXAM MODE

- ① **EXAM:** Select the “NEW” tab located at the top right side of the database display screen.
- ② **Verify Examiner** [Required], **Attending Phys** [Required], and **Referring Phys** (if needed) are entered.
- ③ **ENTER THE EXAM TYPE** [Bio A] and Laterality [OD/OS]. [If the exam is to be used as a teaching case, place a checkmark in the box labeled “Teaching Case”]. Press “OK”
- ④ **ENTER PARAMETERS:** K-readings can be entered for either/both eyes. If the patient has had previous refractive surgery, pressing the “Post Refractive” tab will allow the user to select the formula and IOL Correction [Aramberri Double K, Latkany Myopic/Hyperopic] to apply to the lens power calculation. Then select the “OK” tab.

VERIFY CALIBRATION

- The **Verify Calibration** Window should be displayed if the “**Bio A**” tab had been selected previously.
 - **Select YES** to perform the calibration verification. The VuPad will start emitting an audio “beeping”. Place a drop of water onto the face of the a-scan probe and press the probe onto the calibration cylinder (located on left side probe holder). Once the correct pattern is accepted by the instrument, the beeping will stop and the message “Calibration Verified” will appear.
 - [Note: It may be necessary to increase the “Gain” control in order to achieve the desired pattern].
 - Select NO** to bypass the calibration verification procedure.

EXAM MODE

- Select the desired **OD/OS** tab located on the top of the display. Then proceed as follows: (see Fig. 4)
 - ***SCAN** Tab: Select the desired mode (Cataract, Aphakic, etc.) for the eye to be examined. The user can also select “immersion” mode if using an Immersion Shell for biometry.
 - * **BIOMETRY** Tab: will display the collected data as the scans are performed and accepted.
 - * **SETTINGS** tab: will display the various velocities used for the selected mode. The user can adjust these values if needed by using the ▲ /▼ symbols.




Fig. 4

- The “CONTINUOUS CAPTURE” tab, if selected, will allow the user to set the desired number of scans that the VuPad will obtain with a single press of the footswitch. If “Immersion” mode has been selected, TWO footswitch presses would be required to begin obtaining scans – once to activate the transducer for alignment and one to activate the algorithm for scan acceptance.

- * **HISTORY** tab: records all previous a-scan exams [Date/Time/ Exam type] for the selected patient.

- ② For *“Direct Contact”* scanning, press the *“START SCAN”* tab at the bottom of the display (or press the footswitch). The VuPad will emit a continuous beeping sound indicating that the scan mode has been activated. Place the probe in contact with the cornea and instruct the patient to look at the RED FIXATION light at the tip of the probe. Once the proper alignment is achieved, the VuPad will begin to accept scans which will be shown on the right side of the display screen. Once the desired number of scans are accepted, the VuPad will cease emitting the beeping sound. For *“Immersion”* scanning, the user should press the *“ALIGN PROBE”* tab (or footswitch) until alignment is achieved. Once the user is satisfied that the probe is properly aligned, pressing the *“START SCAN”* tab (or footswitch a second time) will activate the scanning algorithm and begin to collect scan data.
- ③ Selecting the BIOMETRY tab, will permit the user to view all the data collected from the scan session including the ACD, LENS THICKNESS, VITREOUS DEPTH AND AXIAL LENGTH for all scans in the set.
 - The user can omit a scan from being included in the Average by un-checking the *“Enabled”* box located to the left of each data set. This does not delete the scan from the record.
 - The user can *“Delete”* a scan from the record by selecting and holding the individual scan data (for approximately 2 seconds) located on the right side of the display and selecting *“delete”* from the pop-up menu. This will delete the scan and all associated data from the record.

CALCULATION MODE

Once the user has reviewed and is satisfied with the scan results, the CALCULATIONS tab, located at the top of the measurement display screen, should be selected. (See Fig. 4)

- ① Enter the *“Target”* refraction and K-readings (if not entered previously)
- ② Average AXL and ACD should populate based on previous screen.
- ③ **LINK** tab will enable to user to display the same Formula or the same Lens for all 4 of the calculation tables.
- ④ **LENS** can be selected from the drop-down menu which will include all lens models previously added to user profile.
- ⑤ **Formulas** can be selected from the drop down menu which will include all formulas previously added to profile.
- ⑥ IOL table will display with the *“Target”* refraction entered previously at the center. The user can edit the center value by changing the value in the *“Target”* refraction field. (step 1 above).
- ⑦ Once all values are approved and accepted, the user should select the **“ADD TO REPORT”** tab located on the bottom right of the display.
- ⑧ If the fellow eye is also to be measured and included onto the report the user should select either the OD or OS tab at the top of the display and repeat the measurement/calculation steps described above.
- ⑨ Once the eye(s) is measured, calculations performed and **“ADDED TO REPORT”** the **REPORT** tab located at the top of the display should be selected.

REPORTS

Once the REPORT tab is selected, the display will show the report template (see Fig. 5) which includes the A-scan(s) and IOL calculations for each eye for which the measurement/calculations were performed. If required, the user can return to the calculation tables, edit the lens and/or formula data and select the ADD TO REPORT tab again. This will generate a second page to the report which will now display 2 of 2, 3 of 3, etc.....

INCLUDE HEADER: If this tab is highlighted, the default header [SonomedEscalon] will be displayed. Deselecting this tab will cause the Header to be removed from the report.

REMOVE PAGE: Highlighting this tab will remove the currently displayed page from the report.

PRINT: Selecting this tab will display the Windows print window to appear from which the user can select any of the installed printers.

EXPORT: Selecting this tab will display the Export pop-up window to appear for selection of export format by the user.

Press **“CLOSE EXAM”** when finished to return to the patient database screen.

The screenshot shows the SonomedEscalon software interface. At the top, there are two eye sections: OD (Right Eye) and OS (Left Eye). For OD, the measurements are ACD: 2.79 mm, Lens: 5.02 mm, Vitr: 15.65 mm, and AXL: 23.45 mm. For OS, the measurements are ACD: 2.83 mm, Lens: 4.92 mm, Vitr: 15.69 mm, and AXL: 23.45 mm. The interface has a navigation menu on the left with buttons for 'Report', 'History', 'Include Header', 'Remove Page', 'Print', and 'Export'. The main display area shows the 'Report' tab selected, displaying patient information, A-scan graphs, and a grid of IOL calculation tables. The tables include columns for 'Formula', 'IOL', and 'Refraction'. The bottom right corner has a 'Finish Exam' button.

Fig. 5

DIAGNOSTIC A-SCAN MODE

① **Power 'On' System:** By pressing the front panel power switch located on the bottom right side of the unit's front panel.

② **Starting Application:** From Windows desktop, double click on the 'VuMaxHD Software' icon (Fig. 6).




Fig. 6 - VuPad Software Icon

③ **Patient Database Screen:** Upon launching the application the 'Patient Database' screen (Fig. 7) will display. Here you can add a new patient or find previously performed exams.




Fig. 7- Patient Database Screen

④ **Diagnostic A-Utility:** To start the Diagnostic A-Scan application, select the 'Diagnostic A-Utility' button (Fig. 8).




Fig. 8 - 'Diagnostic-A Utility' Button

⑤ **Application Window:** Upon selecting the 'Diagnostic-A Utility' button the VuMAX HD patient database will momentarily close and the Diagnostic A-Scan Utility screen will open (Fig. 9).




Fig. 9 - Diagnostic-A Utility Screen

⑥ **Begin Scan, Stop Scan and Save:** To commence a live scan, select the 'Begin Scan' button. Upon selection an A-Scan trace will appear on the main window. To stop scan, click on the correspondent icon as well as for save a screenshot (Fig. 10)




Fig. 10- 'Begin Scan' 'Stop Scan' and 'Save" Buttons

⑦ **Save As:** Upon selecting the 'Save Screenshot' button, a Windows 'Save As' screen (Fig. 11) will be displayed. Navigate through the folders to select a destination folder for saving the patient's saved A-Scan. Note that a 'File Name' is required for saving the document.




Fig. 11 – 'Save-As' Window

⑧ **Close and Return:** Upon completing examination and saving a record of the test, select the 'Close and Return' button to return to the VuMax HD Patient Database screen.

⑨ **Calibration:** Clicking on the Sonomed Escalon logo while having the control key pressed, it will pop up the S curve calibration as in Fig.12.




Fig. 12 Calibration Screen Window

PACHYMETRY

SELECT EXAM MODE

- ① **EXAM:** Select the “NEW” tab located at the top right side of the database display screen.
- ② **Verify Examiner** [Required], **Attending Phys** [Required], and **Referring Phys** (if needed) are entered.
- ③ **ENTER THE EXAM TYPE** [Pach] and **Laterality** [OD/OS]. [If the exam is to be used as a teaching case, place a checkmark in the box labeled “Teaching Case”]. Press “OK”

PERFORM EXAM




Fig. 13

- ① Select desired exam mode option: (Figures 13 and 14)
 - **Single Point | Single Reading** Provides for a single reading at a single point on cornea.
 - **Single Point | Multiple Readings** Provides for up to 5 readings at a single point on cornea.
 - **Multiple Points | Single Reading** Provides for a single reading at up to 9 points on cornea.
 - **Multiple Points | Multiple Readings** Provides for up to 5 readings at up to 9 points on cornea.Note: If multiple points are to be scanned, select appropriate location on corneal map to first measure.

- ② Confirm Tissue Velocity is appropriate; tap to change as desired.
- ③ Enter measured IOP if corrected IOP based on corneal thickness is desired.
- ④ Select whether a single or multiple points should be automatically taken once system scan criteria are met.
- ⑤ Select “SCAN” button to start scan (or press foot pedal).
- ⑥ Repeat, as necessary, to measure all desired points. Readings can be cleared and repeated, if necessary.

OD
OS

MODE

Single Point | Single Reading

Single Point | Multiple Readings

Multiple Points | Single Reading

Multiple Points | Multiple Readings

TISSUE VELOCITY

1641 m/s

IOP

Measured IOP: Not Entered

IOP Correction: +4 mmHg

Corrected IOP: Not Available

IOP Correction Formula: Ehlers et al
[Click Here to View IOP Correction Table](#) or to
[Select Different IOP Correction Formula](#)

Pachymetry Report

Attending Physician: Escalon, Sonomed

Jones, Michelle

OD

Monday Aug 6 2018 11:50:52.30 AM
Tissue Velocity: 1641 m/s

Measurement Point Results

5 Readings

Average: 486 μ m


Std Dev: 002 μ m

IOP Correction: Not CCT

Point 4 | Reading 4

Measure: 0.485 μ m

Std Dev: 0.003 μ m



Start Scan

READINGS

1
Measure 485 μ m Std Dev 003 μ m
2
Measure 485 μ m Std Dev 001 μ m
3
Measure 490 μ m Std Dev 002 μ m
4
Measure 485 μ m Std Dev 003 μ m
5
Measure 485 μ m Std Dev 002 μ m

Clear Single Reading

Clear All Readings


Perform Probe Test

Finish Exam

Fig 14

REPORT

① Select "Report" button to generate pachymetry report (Fig. 15).



Ideal Vision
123 Main Street
Anytown, USA 12345
215-514-8956

Patient ID	105481	Report Date	Sep 07 2018
Patient Name	Jones, Michelle	Attending Physician	Shatalov, Oleg
Date of Birth	Nov 13 1947	Examiner	Shataloc, Oleg
		Exam Range	Aug 08 2017 – Aug 08 2017


These are general comments regarding the exam...

OD

Monday Aug 6 2018 11:50:52.30 AM
Tissue Velocity: 1641 m/s

Measured IOP 20.1 mmHg
IOP Correction +4 mm Hg
Corrected IOP 24.1 mmHg

Using Ehlers et al IOP Correction Formula




Reading	1	2	3	4	5	Avg	Std Dev	
Point 1	Measure 485	485	490	485	485	486	002	CCT
	Std Dev 002	001	001	002	003			
Point 2	Measure 485	485	490	485	485	486	002	☹
	Std Dev 002	001	001	002	003			
Point 3	Measure 485	485	490	485	485	486	002	
	Std Dev 002	001	001	002	003			
Point 4	Measure 485	485	490	485	485	486	002	☹
	Std Dev 002	001	001	002	003			
Point 5	Measure 485	485	490	485	485	486	002	☹
	Std Dev 002	001	001	002	003			
Point 6	Measure -	-	-	-	-	-	-	Limbus
	Std Dev -	-	-	-	-			
Point 7	Measure -	-	-	-	-	-	-	Limbus
	Std Dev -	-	-	-	-			
Point 8	Measure -	-	-	-	-	-	-	Limbus
	Std Dev -	-	-	-	-			
Point 9	Measure -	-	-	-	-	-	-	Limbus
	Std Dev -	-	-	-	-			

OS

Monday Aug 6 2018 11:50:52.30 AM
Tissue Velocity: 1641 m/s

Measured IOP 20.1 mmHg
IOP Correction +4 mm Hg
Corrected IOP 24.1 mmHg

Using Ehlers et al IOP Correction Formula



Reading	1	2	3	4	5	Avg	Std Dev	
Point 1	Measure 485	485	490	485	485	486	002	CCT
	Std Dev 002	001	001	002	003			
Point 2	Measure 485	485	490	485	485	486	002	☹
	Std Dev 002	001	001	002	003			
Point 3	Measure 485	485	490	485	485	486	002	
	Std Dev 002	001	001	002	003			
Point 4	Measure 485	485	490	485	485	486	002	☹
	Std Dev 002	001	001	002	003			
Point 5	Measure 485	485	490	485	485	486	002	☹
	Std Dev 002	001	001	002	003			
Point 6	Measure 472	470	471	471	471	471	001	Limbus
	Std Dev 001	002	001	001	002			
Point 7	Measure 468	465	462	462	468	465	003	Limbus
	Std Dev 002	002	003	003	002			
Point 8	Measure 472	470	471	471	471	471	001	Limbus
	Std Dev 001	002	001	001	002			
Point 9	Measure 468	465	462	462	468	465	003	Limbus
	Std Dev 002	002	003	003	002			

1/1

Fig 15

19

B-SCAN AND UBM

CONFIGURE SYSTEM

- ① After logging in by entering the administrator name and password, the home screen is the patient database, which shows a history of patients and their respective exams. From this screen, select the <Configure> button.
- ② **ENTER ASSOCIATE INFORMATION** by selecting the <New> button, entering information and clicking the <Apply> button. Each associate can be designated as an examiner, attending physician, and/or referring physician.
- ③ **ENTER ULTRASOUND INFORMATION** by selecting the Ultrasound tab and selecting the <Add...> button. Enter the frequency in MHz of any water path transducer(s).
- ④ **ENTER DIAGNOSES FOR REFERENCE ON EXAM RECORDS** by selecting the Diagnoses tab and selecting the <New Diagnosis> tab.

PATIENT DATABASE









Upon opening the VuPad program, the home screen shown is the patient database, which shows a history of patients and their respective exams. From this screen, several functions can be performed:

- ① **SEARCH** for patients by entering the known information and selecting the <search> button. A list of patients matching the criteria will be shown. Select the <clear> button to return to a list of all patients.
- ② **ENTER NEW PATIENT** by selecting the <new patient> button, entering appropriate information, and selecting the <Apply> button.
- ③ **DELETE PATIENT** by selecting the <Delete> button in the patient section.
NOTE: Deleting a patient record will also delete all associated exam records. A prompt will appear to confirm intent to actually delete a patient record.
- ④ **SELECT PATIENT** by highlighting the patient from the patient list – a history of any saved exams for the patient will appear in the exam list.
- ⑤ **SELECT EXISTING EXAM** by double-clicking on an exam record for a particular patient. This function opens the exam and displays associated images (see EXAM RECORD section of this user manual).
- ⑥ **CREATE NEW EXAM** by selecting the <New Exam> button.
NOTE: Patient information must be entered first.
- ⑦ **DELETE EXAM RECORD** by selecting the <Delete> button in the exam section.
NOTE: Deleting an exam record will delete all associated images and data. A prompt will appear to confirm intent to actually delete an exam record.
- ⑧ **IMPORT EXAM RECORDS** by selecting the <Import> button and following the on-screen instructions.
- ⑨ **EXPORT EXAM RECORDS** by selecting the <Export> button and following the on-screen instructions.
NOTE: This function may be useful when connecting to an EMR system, referrals, or creating presentations.
- ⑩ **PRINT REPORT** by selecting the <Report> button. In the Report screen, the images may be selected (across the patient's history of exams) and the format chosen for one, two, four, or six images per page.
NOTE: The report header may be toggled on or off by selecting the <Header> button.
- ⑪ **CONFIGURE SYSTEM** by selecting the <Configure> button to enter examiner, physician, and referring physician information, among other set-up tasks.

EXAM MODE

After selecting the <New Exam> button from the Patient Database screen, enter the exam information, the Exam screen is shown from which imaging can be performed.

- ① **SELECT THE PROBE FREQUENCY** to be used from the drop-down menu and select the <Refresh Probes> button.
 Standard B-Scan Mode: 12 MHz or 20 MHz probes may be used.
 UBM Mode: 35 MHz or 50 MHz transducer may be used in conjunction with a water path probe.
- ② **SELECT APPROPRIATE PRESET SCAN MODE** and/or adjust image controls for optimizing image quality for the particular structure of interest. Available image controls and preset scan modes are shown on the next page.
NOTE: Image controls may also be adjusted during or after a scan.
- ③ **PLACE PROBE INTO POSITION AND START SCANNING** by either depressing the footswitch or selecting the <New Scan> button. When the footswitch is pressed or <Stop Scan> button selected, a video clip of the prior 50 frames is captured.
DANGER! WHEN USING THE WATER PATH PROBE IN UBM MODE (35 OR 50 MHZ), THE USER MUST TAKE EVERY PRECAUTION TO PREVENT THE TRANSDUCER FROM TOUCHING THE EYE. CONTACT BETWEEN A MOVING TRANSDUCER AND THE EYE CAN CAUSE SEVERE INJURY.
- ④ **REVIEW SCAN** by selecting the Play button and selecting the playback rate (or scroll through video clip one frame at a time). Image controls may be adjusted to see the effect on the displayed image, including use of the zoom feature. Individual frames may be saved by selecting the <Save Frame> button.
- ⑤ **ANNOTATE SCAN** by selecting Tools tab and selecting the desired function. Annotations are saved with the associated video clip or image.
NOTE: A frame is automatically saved when an image is annotated.

ANNOTATIONS	
	Angle Measurement places adjustable angle edges onto image and displays angle in degrees.
	Arbitrary A-Scan allows user to see an A-scan trace along any vector.
	Linear Measurement places adjustable calipers onto image and displays distance.
	Pan allows user to move image while holding down the left mouse button.
	A-Scan displays an A-scan trace along the axis of a B-scan image. (B-scan only)
	Measurement Template displays approximate common measurements for UBM sulcus-to-sulcus video clips with ability to move measurement endpoints. (UBM-only)
	Refresh Measurements will reset measurements to those auto-generated by system. (UBM-only)
	Clear removes the highlighted annotation from the image.
Diagnosis	Add or Remove a diagnosis by selecting the appropriate buttons.
Notes	Add or Remove notes by selecting the appropriate buttons.

- ⑥ **EXPORT FRAME** by selecting Tools tab and selecting <Export JPG> button.
- ⑦ **USE NOMOGRAMS TO DETERMINE ICL LENS SIZING** by selecting the Tools tab and selecting the desired nomogram and entering required information to yield the resultant lens size.
- ⑧ **EXIT EXAM MODE** by selecting <Finish Exam> button

IMAGE CONTROLS	
TVG	Time Variable Gain (TVG) attenuates near field portion of scan sector.
Baseline	Affects low-level signal threshold to help reduce noise.
Zoom	Adjusts area of scan sector display on monitor.
Log Gain	Adjusts overall amplification of received signal.
Exp Gain	Adjusts overall signal amplitude in image with minimal effect on low-level noise (otherwise known as “E-Gain”).

UBM PRESET MODES	
Sulcus-to-Sulcus	Optimized to measure sulcus-to-sulcus for lens sizing with wide 30° scan angle. Note: When measuring sulcus-to-sulcus, ensure iris plane is saturated, with anterior and posterior lens echoes visible. It is advisable to perform multiple scans with measurements to ensure most accurate results.
Angle Detail	With scan angle of 15° and mid-range settings for gain and scan frequency, optimized to provide excellent detail of angle structures, tumors, traumas, and fine details.

B-SCAN PRESET MODES	
Orbit	Optimized to image the orbit with 50° scan angle.
Deep Retina / Choroid	Optimized to view beyond the surface of the retina and into the choroid, with a 50° scan angle.

MAINTENANCE *No preventive maintenance required.*

Always inspect probe tips prior to use to ensure no scratches or damage are present. Do not use probes if damage is visible.

SYSTEM CLEANING

Clean the VuPad™ probes with a damp cloth to remove gel and other debris. Use appropriate products to clean the cables and foot pedals as necessary.

WARNING: Disconnect the AC POWER before cleaning the system.

PROBE CLEANING AND DISINFECTION

The probe must be cleaned and disinfected between patients to prevent patient-to-patient transmission of infection. Prior to any cleaning or disinfecting, unplug the probe from its cable. Cleaning is intended to remove dirt and debris from the probe, and to reduce the presence of microorganisms. Disinfection is performed after cleaning to address microorganisms.

PROBE CLEANING

- a) A few drops of common concentrated dishwashing detergent or enzymatic product diluted in warm tap water may be used. Scrub the probe in a soapy solution that facilitates the suspension and washing away of the unwanted contaminants. The probes may be vigorously scrubbed, as needed, to remove contaminants. A soft bristle brush may be used to scrub the narrow gap where the probe window joins the probe cover.
- b) Rinse the probe thoroughly with distilled or deionized water and allow to air dry or blot dry with a clean, soft, lint-free cloth that does not leave lint or debris
- c) If a water path probe is used, the transducer may be washed with light pressure. Do not scrub the face of the transducer harshly or use any abrasive cleaners or cloths. Careful cleaning will lengthen the useful life of the transducer. *Special care should be taken not to rub the gold surface of the transducer.*

PROBE DISINFECTION

- a) For low-to-moderate disinfection: After cleaning, immerse the probe in 70% isopropyl alcohol (70% IPA) for 5-10 minutes. Rinse the probe thoroughly with distilled or deionized water and allow to air dry or blot with a clean, soft, lint-free cloth that does not leave lint or debris. If not used immediately, the probe may be placed in a clean bag for storage.
- b) For a higher level of disinfection: After cleaning, immerse the probe in 2-3% W/W hydrogen peroxide for 8-10 minutes. Rinse the probe thoroughly with distilled or deionized water and allow to air dry or blot with a soft, lint-free cloth or gauze. If not used immediately, the probe may be placed in a clean bag for storage.
- c) Probes are hermetically sealed and, if necessary, the entire probe (up to the connector) may be immersed in disinfecting solution. However, this should be reserved for rare cases where it is judged by the clinician that the entire probe has been contaminated. When this is not the case, then it is only necessary to immerse the portion of the probe that has been in contact with the patient, plus approximately 2 cm. of the probe cover.

CAUTIONS: NEVER IMMERSE THE CONNECTOR AT THE END OF THE PROBE!

NEVER AUTOCLAVE ANY PROBES OR TRANSDUCERS OR EXPOSE TO HIGH HEAT; EXCESSIVE TEMPERATURE WILL CAUSE DAMAGE

PROBE CLEANING & DISINFECTION RECOMMENDATIONS

- Do not allow probes to come in contact with any solutions for longer than 10 minutes at a time. A longer contact time is at the discretion of the user.
- Thorough rinsing with distilled or deionized water is recommended after contact with any cleaning or disinfectant agent to remove traces of the solution.
- Air drying is acceptable following water rinse.
- If a cloth is desired for drying, blot dry with a clean, soft, lint-free material that does not leave visible debris or lint on the probe.
- Diluted sodium hypochlorite solution (1:10 bleach solution) may be used as a high-level disinfectant on probes that do not have an external transducer. Always rinse very thoroughly and do not use for longer than 10 minutes.
- Do not use bleach on UBM transducers. The gold plated surfaces of the UBM transducers should not come in contact with bleach.
- FDA cleared disinfectants for low-level or high-level disinfection may be used on the probes according to the facility and/or manufacturer's instructions, but not longer than 10 minutes. A longer contact time is at the discretion of the user.
- Thorough and continuous rinsing with copious amounts of the disinfectant solution for several minutes is an effective means of disinfecting between patients.
- Avoid use of abrasives on all probes.

CAUTIONS WHEN CLEANING WATER-PATH PROBES

- If necessary, the transducer may be detached from the water path probe and the two disinfected separately. If detached for cleaning, take care that the surface mount connector (SMC) of the probe and transducer are thoroughly dried before remounting the transducer to the probe. The internal surfaces of the transducer connections must be free of moisture prior to re-assembly. *Special care should be taken not to rub the gold surface of the transducer.*
- The hard, plastic parts of the probe may be scrubbed, as needed to remove contaminants. A soft bristle brush may be used to scrub crevices. Take care to thoroughly, but only lightly, scrub the soft rubber part surrounding the transducer connector.
- The gold surface of the transducer may be washed with very light pressure. Take care NOT to scrub the face of the transducer harshly or to use any abrasive cleaners or cloths. Careful cleaning will lengthen the useful life of the transducer.
- *Sodium Hypochlorite* (i.e. Bleach): Diluted (0.6% concentrations) and undiluted bleach solutions attack the gold plated surfaces of the transducer, and once a breach has been made, the sensitivity of the transducer may be diminished (proportional to the loss of gold plating). Therefore, methods employing 0.6% concentrations, or more, of sodium hypochlorite should be avoided whenever possible.
- *Tap water*: Contaminants in common tap water may contribute to corrosion of the metal surfaces of the probe and transducer after extended use. The use of tap water, rather than deionized water or distilled water is, therefore, contraindicated. If tap water is used, rinse thoroughly with distilled or deionized water and lightly wipe the surface clean of any residues.
- *Abrasives*: The gold plating on the transducer surface is very thin and cannot be thickened without negatively affecting the transducer frequency and performance. Therefore, extreme care must be taken to avoid scratching or otherwise damaging the gold surface of the film. Never use any kind of abrasive cloth or tissue when wiping the transducer surface. Camera lens cleaning paper or soft gauze may be used. Apply no more force than is needed when wiping the transducer surface. Minor scratches to the gold surface will not damage performance. Contact the manufacturer if the gold plating is breached and a silver color is seen in place of the gold.

Scan Type	Probe Type	Use	Patient Contact	Device Classification	Disinfection Level	Recommended	Method
A Scan / Pachymetry	A Probe / Pachymeter Probe	Direct Contact	Mucous Membrane	Semi-Critical	High Level	2-3% Hydrogen Peroxide	Clean probe as needed. Soak for 8-10 minutes. Rinse thoroughly with distilled or DI water.
A Scan	A Probe	Immersion	None	* Semi-Critical	High Level	2-3% Hydrogen Peroxide	
B Scan	B probe	Over the Eyelid	Intact Skin	Non-Critical	Low Level	70% Isopropyl Alcohol	Clean probe as needed. Soak for 5-10 minutes. Rinse thoroughly with distilled or DI water.
B Scan	B probe	Direct Contact	Mucous Membrane	Semi-Critical	High Level	2-3% Hydrogen Peroxide	Clean probe as needed. Soak for 8-10 minutes. Rinse thoroughly with distilled or DI water.
UBM	Water Path Probe with External Transducer	Clear Scan Cover	None	NA	NA	Cleaning or low level disinfection if necessary	Clean probe as needed. 70% Isopropyl alcohol wiped on the probe. Rinse thoroughly with distilled or DI water.
UBM	Water Path Probe with External Transducer	Immersion	None	* Semi-Critical	High Level	2-3% Hydrogen Peroxide	Clean probe as needed. Soak for 8-10 minutes. Rinse thoroughly with distilled or DI water.

* The immersion technique does not involve patient contact, however, the mucous membrane is exposed to the immersion solution and the probe is used in the immersion solution. Therefore, worst case conditions are considered in assigning the device classification as semi-critical.
DI Water = Deionized Water

EYE CUP AND PRAGER SHELL CLEANING AND DISINFECTION

- Eye cups are manufactured using acrylic and silicone. The eye cups may be cleaned using the same procedure as the probes using water with detergents or enzymatic products. Cleaning is performed as needed to remove foreign material (e.g., soil, and organic material).
- The use of facility disinfectants may be used to disinfect eye cups according to label instruction and/or facility procedures followed by thorough distilled or DI water rinse and air dry. Sonomed does not prescribe specific instruction for disinfection of eye cups, however, the eye cups should not be subjected to autoclaving or extreme heat.
- Follow the instructions provided by ESI, Inc. included with the Prager Shell. The Prager tray may be used for cleaning and disinfecting. Place the shell into the tray with the Luer fitting and setscrew facing the notched area of the tray. If used as a cleaning tray for disinfection, fill tray approximately 2/3 with an antiseptic solution (ESI references 7.5 % H₂O₂). After removing the shell from the tray, thoroughly rinse and let air dry.

DISPOSAL

When disposal is required, the equipment and associated cleaning and disinfecting chemicals should be disposed of in accordance with local, state, and federal laws.

In the European Union, follow Waste Electrical & Electronic Equipment (WEEE) Directive 2012/19/EU Annex I, 4.07.2018

TECHNICAL SERVICE AND SUPPORT

There are no user-serviceable parts within the system. Please contact Sonomed Escalon or your local distributor to request technical service and support. Technical support 800-227-1285 or 516-354-0900.

Email: ultrasound-support@escalonmed.com

EVENT REPORTING

In the event a serious injury or incident occurs in relation to the use of the ophthalmic ultrasound, the event should be reported to Sonomed Escalon and the applicable regulatory authorities.

SPECIFICATIONS / PERFORMANCE CHARACTERISTICS

VuPad Unit General Specifications	
Operating System	Windows 10 Enterprise
Processor	Pentium N4200, 1.1 GHz (2.0 GHz Turbo) Quad Core
Memory / Date Storage	8GB / 500 GB SSD (1 TB optional)
Connectivity	2 x USB 3.0, Ethernet, HDMI, WiFi, Bluetooth 4.0
Dimensions	13.3" W x 8.0" L x 2.0" D (33.8 cm x 20.32 cm x 5.08 cm)
Weight	4.5 lbs (2.1 kg)
Power Supply	Input: 100 – 240VAC, 50/60 Hz Output: 12VDC / 5A
Device	Touchscreen (optional wireless keyboard / mouse)
Display Screen	Built-in, 10.1" High Resolution Touchscreen LCD with optional simultaneous display on external LCD monitor via HDMI port
Display Resolution	1280 x 800 HD
Display Mounting	Adjustable kickstand (optional articulate monitor arm)
Printer	Any Windows 10 compatible printer, plus Sony/Mitsubishi video printer
Operating Conditions	Operating Temperature: 41 to 104°F (5 to 40°C) Storage Temperature: -40 to 158°F (-40 to 70°C) Operating Relative Humidity: 0 to 90% non-condensing Storage Relative Humidity: 0 to 90% non-condensing Atmospheric Pressure: 70 to 101kPa
Wireless Specifications	Frequency band of operation: 5 GHz Operating Frequency: 5150 – 5350 MHz Bandwidth: 802.11 a/n/ac: 20/40/80/MHz RF Output Power: 23 dBm max. (5150 – 5725 MHz) IEEE 802.11 a/n/ac

A-SCAN SPECIFICATIONS	
Transducer Frequency	10 MHz
Active Diameter	3.5 mm
Examination Modes	Cataract, Dense Cataract, Phakic, Aphakic, Pseudophakic, Manual Mode
Gate Adjustment	Automatic and manual 4-gate positions with auto-detection of scleral spike
Adjustable GAIN Range	40 – 80 dB
Electronic Resolution	± .023 mm
Clinical Resolution	± 0.10 mm
Axial Resolution	0.019 mm
Scan Depth	45 mm

A-SCAN SPECIFICATIONS (cont'd)	
Focal Length	25 mm
Measurement Range	18 – 40 mm
Auto-Calibration	YES
Method	Contact and Immersion
Capture Mode	Automatic, Manual
Tissue Velocity	Adjustable velocity for each tissue segment
Biometry Measurement Parameters	Anterior Depth (ACD) Lens Thickness (LT) Vitreous (VD) Axial Length (AXL) Automatic calculation of average axial length and standard deviation based on 10 measurements
IOL Formulas	Theoretic – T (SRK-T), Regression II (SRK-II), Hoffer-Q, Binkhorst, Holladay, Haigis
Post-Op Refractive Formulas	Latkany Myopic Regression, Latkany Hyperopic, Aramberri Double-K
Diagnostic A-scan	Optional Diagnostic A-scan module with 8 MHz diagnostic A-scan probe

B-SCAN SPECIFICATIONS	
Transducer Frequency	12 MHz (<i>optional 20 MHz</i>)
Scan Angle	50 degrees (<i>40 degrees</i>)
Scan Depth	60 mm
Focal Zone	13 – 35 mm
Focal Point	24 mm (<i>22 mm</i>)
Axial Resolution	0.013 mm (<i>0.0095 mm</i>)
Lateral Resolution	0.068 mm (<i>0.060 mm</i>)
Frame Rate	20 fps
Grey Levels	256
Time Variable Gain	0 – 7 dB
Adjustable GAIN Range	40 – 130 dB
Exponential Gain (EGain)	Enhanced gain adjustment of 0 – 100 dB
Examination Presets	Deep Retina/Choroid, Vitreous Body, Retina Surface, Orbit and Silicon Oil
Measurements & Annotations	Caliper, Arbitrary A-scan, Angle, Area, Text Box, Marker

PACHYMETER SPECIFICATIONS	
Transducer Frequency	20 MHz
Measurement Range	0.3 mm – 1 mm
Clinical Accuracy	± .005 mm
Electronic Accuracy	± .001 mm
Sample Rate	192 MHz
Bias Correction	80 – 120 %
Tissue Velocity	Adjustable
Method	Contact
Measurements	Auto-Capture, Continuous and Sensing Algorithm Auto Calibration and Probe Test CCT and cartographic map
Corneal Map	Multiple graphical corneal maps: Selectable central corneal thickness and various peripheral areas
Scan Modes	Single Point – Single Reading Single Point – Multiple Readings Multiple Points – Single Reading Multiple Points – Multiple Readings
IOP Correction	CCT-based automatic IOP correction using formulas: Ehlers et al., Doughty & Zaman and Pillunat et al.

UBM SPECIFICATIONS	
Transducer Frequency	35 MHz (<i>optional 50 MHz</i>)
Scanning Angle	30 degrees (max)
Frame Rate	20 fps (max)
Max Scan Depth	22 mm
Focal Depth	12 mm
Focal Zone	9 – 13 mm
Axial Resolution	.007 mm (<i>.004 mm</i>)
Lateral Resolution	.039 mm (<i>.035 mm</i>)
Adjustable GAIN Range	40 – 130 dB
Exponential Gain (EGain)	Enhanced gain adjustment of 0 – 100 dB
Time Variable Gain	0 – 7 dB
Examination Presets	Sulcus-to-Sulcus, High Resolution, Angle Analysis, Motion Picture

DATA AND IMAGING PROCESSING	
Number Saved Images (Frames)	Unlimited
Video Clip Length	12 video clips per exam (50 frames per video clip)
Video Clip Replay	Real-time, scalable slow-motion, frame-by-frame
Number Saved Video Clips per Eye	6 per exam (unlimited exams)
Default Annotation (Frame)	Patient name, ID, user ID, exam date, DOB, Position scan indicator, Gain, E Gain, baseline, eye tracking overlay, probe frequency, scan angle
Additional User Annotation (Frame)	Caliper and angle dimensions, user defined free-form text, arbitrary A-scan
Export File Format	.jpg (images) and .avi (video clips)
Report Template	User customizable report template
Export Report Format	.pdf and DICOM
Connectivity	802.11n dual-band WiFi, Bluetooth 4.0, GigE Ethernet LAN, USB 3.0 ports DICOM compatible (optional license)
Image Brightness / Contrast	Adjustable: TVG, Log, Baseline and Gain
Zoom Magnification	Continuous Interpolative Zoom (0.5x – 4.0x)
Measurement Calipers	Unlimited: Line/Distance, Angle, Area and Arbitrary A-scan Caliper
Eye Tracking Algorithm	Automatic alignment tracking with measurements of CCT, ACD, LT, ATA, STS
Quantitative Angle Analysis Software	Complete and fully automatic angle analysis for AOD, TIA, ARA, TISA

ALARA SECTION AND EMISSIONS

("As Low As Reasonably Achievable")

STATISTICAL ANALYSIS OF MEASURED DATA

A statistical analysis was performed on the data to examine the upper output limits based on a one-sided tolerance limit approach. The mean and standard deviation of the Spatial-Peak, Time-Average Intensity and Mechanical Index were found, and the upper output limits were calculated from the following formula:

$$X = x + K * Sx$$

Where X is the upper output parameter limit, x is the average of the measured output parameter, and Sx is the standard deviation of the measured output parameter. A value of K was chosen which corresponds to a 90% probability that 90% of all probes would fall below the calculated limits of X.

RESULTS

Statistical analysis showed that the probes tested produced MI and Ispta.3 values below FDA limit values.

ACCURACY

The accuracy of the emissions figures is approximately 26.6% for all intensity values reported, 13.3% for all pressure values reported and 13.3% for the Mechanical Index.

CAUTION

Care should be taken by user to minimize exposure of patient to ultrasound energy by keeping examinations as short as possible.

Probe	35 MHz Transducer		50 MHz Transducer	
Mode	B-Mode (UBM)		B-Mode (UBM)	
Material	Gold-plated polymer membrane		Gold-plated polymer membrane	
Nominal Center Frequency	35 MHz		50 MHz	
Pulse Repetition Frequency	2560 Hz		2560 Hz	
Type	B-Scan (energy emitted during active scan)		B-Scan (energy emitted during active scan)	
Measure	MI [no units]	$I_{SPTA,3}$ [mW/cm ²]	MI [no units]	$I_{SPTA,3}$ [mW/cm ²]
Sample Size	3.0	3.0	3.0	3.0
K	4.258	4.258	4.258	4.258
Mean	0.086	0.025	0.023	0.002
Standard Deviation	0.0070	0.0030	0.0040	0.0004
Limit	0.115	0.037	0.038	0.003

Probe	12 MHz Transducer		20 MHz Transducer	
Mode	B-Mode		B-Mode	
Material	TPX Polymer		TPX Polymer	
Nominal Center Frequency	12 MHz		20 MHz	
Pulse Repetition Frequency	2560 Hz		2560 Hz	
Type	B-Scan (energy emitted during active scan)		B-Scan (energy emitted during active scan)	
Measure	MI [no units]	$I_{SPTA,3}$ [mW/cm ²]	MI [no units]	$I_{SPTA,3}$ [mW/cm ²]
Sample Size	3.0	3.0	3.0	3.0
K	4.258	4.258	4.258	4.258
Mean	0.149	0.105	0.135	0.240
Standard Deviation	0.0154	0.0260	0.0117	0.0504
Limit	0.215	0.216	0.184	0.455

Probe	10 MHz (A-Scan) Transducer	
Mode	A-Mode	
Material	Lead Metaniobate	
Nominal Center Frequency	10 MHz	
Pule Repetition Frequency	5880 Hz	
Type	A-Scan (energy emitted during active scan)	
Measure	MI [no units]	I _{SPTA,3} [mW/cm ²]
Sample Size	3.0	3.0
K	4.258	4.258
Mean	0.186	7.82
Standard Deviation	.00458	0.570
Limit	0.206	10.3

Probe	20 MHz (Pachymeter) Transducer	
Mode	A-Mode	
Material	Ceramic	
Nominal Center Frequency	20 MHz	
Pule Repetition Frequency	5880 Hz	
Type	A-Scan (energy emitted during active scan)	
Measure	MI [no units]	I _{SPTA,3} [mW/cm ²]
Sample Size	3.0	3.0
K	4.258	4.258
Mean	0.0683	4.02
Standard Deviation	0.00681	0.809
Limit	0.0973	7.47

NOTE: The energy will always be attenuated by the tissue between the transducer and the focus when used as recommended. The values presented here are the values at the focal point, the point of maximum intensity.

ACOUSTIC OUTPUT REPORTING TABLES FOR TRACK 1: AUTO SCANNING MODE

Note: Track 1 Reporting Tables show the worst-case indices for each probe type and operating conditions that must be reported.

Transducer Model		Sonomed Escalon 10 MHz Probe (s/n D09A655)			
Operating Mode		A-Mode			
Application		Ophthalmic			
Acoustic Output		MI [no units]	I _{SPTA.3} [mW/cm ²]	I _{SPPA.3} [W/cm ²]	
Global Maximum Value		0.190	8.42	8.52	
Associated Acoustic Parameters	Pr.3 [Mpa]	0.531			
	W ₀ [mW]		0.538	0.538	
	f _c [MHz]	7.80	7.80	7.80	
	Z _{sp} [cm]	1.90	1.90	1.90	
	Beam Dimensions	x ₋₆ [cm]		0.177	0.177
		y ₋₆ [cm]		0.162	0.162
	PD [μS]	0.168		0.168	
	PRF [Hz]	5880		5880	
	EDS	Az [cm]		0.47	
Ele. [cm]			0.47		

Transducer Model		Sonomed Escalon 12 MHz Probe (s/n 22064)			
Operating Mode		B-Mode			
Application		Ophthalmic			
Acoustic Output		MI [no units]	I _{SPTA.3} [mW/cm ²]	I _{SPPA.3} [W/cm ²]	
Global Maximum Value		0.167	0.135	12.9	
Associated Acoustic Parameters	Pr.3 [Mpa]				
	W ₀ [mW]		0.0832	0.0832	
	f _c [MHz]	12.2	12.2	12.2	
	Z _{sp} [cm]	1.60		1.60	
	Beam Dimensions	x ₋₆ [cm]			0.0689
		y ₋₆ [cm]			0.0824
	PD [μS]	0.131		0.131	
	PRF [Hz]	2560		2560	
	EDS	Az [cm]		0.600	
Ele. [cm]			0.600		
Operating Control Conditions		Scan angle 60°	256 lines per frame	Scan rate 10 Hz	

Transducer Model		Sonomed Escalon 20 MHz Probe (s/n 22160)			
Operating Mode		B-Mode			
Application		Ophthalmic			
Acoustic Output		MI [no units]	I _{SPTA.3} [mW/cm ²]	I _{SPPA.3} [W/cm ²]	
Global Maximum Value		0.122	0.184	11.5	
Associated Acoustic Parameters	Pr.3 [Mpa]	0.456			
	W ₀ [mW]		0.0296	0.003	
	f _c [MHz]	14.0	14.0	14.0	
	Z _{sp} [cm]	1.70		1.70	
	Beam Dimensions	x ₋₆ [cm]			0.0601
		y ₋₆ [cm]			0.0605
	PD [μS]	0.0650		0.0650	
	PRF [Hz]	2560		2560	
	EDS	Az [cm]		0.7	
Ele. [cm]			0.7		
Operating Control Conditions		Scan angle 60°	256 lines per frame	Scan rate 10 Hz	

Transducer Model		Sonomed Escalon 35 MHz Transducer (s/n 35-00219)			
Operating Mode		B-Mode			
Application		Ophthalmic			
Acoustic Output		MI [no units]	I _{SPTA.3} [mW/cm ²]	I _{SPPA.3} [W/cm ²]	
Global Maximum Value		0.094	0.028	4.780	
Associated Acoustic Parameters	Pr.3 [Mpa]	0.460			
	W ₀ [mW]		0.003	0.003	
	f _c [MHz]	24.095	24.095	24.095	
	Z _{sp} [cm]	1.000		1.000	
	Beam Dimensions	x-6 [cm]			0.039
		y-6 [cm]			0.038
	PD [μS]	0.043		0.043	
	PRF [Hz]	2560		2560	
	EDS	Az [cm]		0.7	
Ele. [cm]			0.7		
Operating Control Conditions		Scan angle 15°	256 lines per frame	Scan rate 10 MHz	

Transducer Model		Sonomed Escalon 50 MHz Transducer (s/n 50-00376)			
Operating Mode		B-Mode			
Application		Ophthalmic			
Acoustic Output		MI [no units]	I _{SPTA.3} [mW/cm ²]	I _{SPPA.3} [W/cm ²]	
Global Maximum Value		0.019	0.001	0.190	
Associated Acoustic Parameters	Pr.3 [Mpa]	0.112			
	W ₀ [mW]		0.0003	0.0003	
	f _c [MHz]	34.168	34.168	34.168	
	Z _{sp} [cm]	1.000		1.000	
	Beam Dimensions	x-6 [cm]			0.035
		y-6 [cm]			0.041
	PD [μS]	0.056		0.056	
	PRF [Hz]	2560		2560	
	EDS	Az [cm]		0.7	
Ele. [cm]			0.7		
Operating Control Conditions		Scan angle 15°	256 lines per frame	Scan rate 10 MHz	

Transducer Model		Sonomed Escalon 20 MHz Pachymeter Probe (s/n PD20.0-0108-0341)			
Operating Mode		A-Mode			
Application		Ophthalmic			
Acoustic Output		MI [no units]	I _{SPTA.3} [mW/cm ²]	I _{SPPA.3} [W/cm ²]	
Global Maximum Value		0.0760	4.91	3.69	
Associated Acoustic Parameters	Pr.3 [Mpa]	0.323			
	W ₀ [mW]		0.236	0.236	
	f _c [MHz]	17.9	17.9	17.9	
	Z _{sp} [cm]	1.00	1.00	1.00	
	Beam Dimensions	x-6 [cm]		0.153	0.153
		y-6 [cm]		0.142	0.142
	PD [μS]	0.0780		0.0780	
	PRF [Hz]	1.70E+4		1.70E+4	
	EDS	Az [cm]		0.2	
Ele. [cm]			0.2		

IMMUNITY TEST LEVELS

Electrostatic Discharge

Air Discharge Voltage: 2.0 kV, 4.0 kV, 8.0 kV, 15.0 kV
 Contact Discharge Voltage: 2.0 kV, 4.0 kV, 8.0 kV

Radiated Immunity, 80 MHz to 2.7 GHz and Proximity Fields

Radiated Immunity, Test Parameters		
Frequency Range	80 to 1000 MHz	1.0 to 2.7 GHz
Field Strength	10 V/m	10 V/m
Window Tested	NA	1
Modulation	1 kHz, 80%, AM	
Dwell Time	0.5 seconds	
Polarization of Applied Field	Horizontal and Vertical	

Radiated Immunity, Proximity Field Parameter			
Frequency (Hz)	Modulation	Level V/m	Test Distance
385	Pulse, 18 Hz, 50% DC	27	1.0 m
450	FM, 1 kHz Sine, ± 5 Hz Deviation	28	1.0 m
710, 745, 780	Pulse, 217 Hz, 50% DC	9	1.0 m
810, 870, 930	Pulse, 18 Hz, 50% DC	28	1.0 m
1720, 1845, 1970	Pulse, 217 Hz, 50% DC	28	1.0 m
2450	Pulse, 217 Hz, 50% DC	28	1.0 m
5240, 5500, 5785	Pulse, 217 Hz, 50% DC	9	1.0 m

Electrical Fast Transient / Burst

Transient Voltage: 0.5 kV, 1.0 kV, 2.0 kV
 Transient Polarity: Positive and Negative
 Repetition Rate: 100 kHz
 Rise Time of Pulse: 5 nS ± 30%
 Pulse Duration: 50 nS ± 30%
 Burst Period: 300 ms ± 20%
 Burst Duration: 15 ms ± 20%

Surge, Power Ports

Voltage: 0.5 kV, 1.0 kV Differential Mode
 0.5 kV, 1.0 kV Line to Line Mode
 0.5 kV, 1.0 kV, 2.0 kV Common Mode
 Polarity: Positive and Negative
 Pulse Phase: 0°, 90°, 180°, 270°
 Open Circuit: Rise Time 1.2 µsec, Duration 50.0 µsec, Rep Rate 1 ppm
 Short Circuit: Rise Time 8.0 µsec, Duration 20.0 µsec

Conducted Immunity Including ISM and Amateur Radio Bands

AC Mains and I/O Cables

Frequency Range, Test Level: 0.15 to 80 MHz, 10 Vms
 ISM Test Frequencies: 6.765 to 6.795, 13.553 to 13.567
 26.957 to 27.283, 40.660 to 40.700 MHz
 Test Level: 6 Vrms

Power Frequency, Magnetic Immunity

Frequency: 50 Hz and 60 Hz
 Applied Signal Level: 30 A/M RMS

VOLTAGE DIPS AND INTERRUPTS

100 VAC 60 Hz							
Specification	Rated Voltage	Frequency (Hz)	Voltage Test Level (%)	Voltage Dip (%)	Test Voltage (Vac)	Duration (Periods)	Result
0% UT for 0.5 Cycles @ 0°	100	60	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 45°	100	60	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 90°	100	60	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 135°	100	60	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 180°	100	60	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 225°	100	60	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 270°	100	60	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 315°	100	60	0	100	0	0.5 Cycle	(1)
0% UT for 1 Cycle @ 0°	100	60	0	100	0	1 Cycle	(1)
70% UT for 30 Cycles @ 0°	100	60	70	30	70	30 Cycles	(1)
0% UT for 300 Cycles @ 0°	100	60	0	100	0	300 Cycles	(2)

220 VAC 60 Hz							
Specification	Rated Voltage	Frequency (Hz)	Voltage Test Level (%)	Voltage Dip (%)	Test Voltage (Vac)	Duration (Periods)	Result
0% UT for 0.5 Cycles @ 0°	220	60	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 45°	220	60	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 90°	220	60	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 135°	220	60	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 180°	220	60	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 225°	220	60	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 270°	220	60	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 315°	220	60	0	100	0	0.5 Cycle	(1)
0% UT for 1 Cycle @ 0°	220	60	0	100	0	1 Cycle	(1)
70% UT for 30 Cycles @ 0°	220	60	70	30	154	30 Cycles	(1)
0% UT for 300 Cycles @ 0°	220	60	0	100	0	300 Cycles	(2)

240 VAC 50 Hz							
Specification	Rated Voltage	Frequency (Hz)	Voltage Test Level (%)	Voltage Dip (%)	Test Voltage (Vac)	Duration (Periods)	Result
0% UT for 0.5 Cycles @ 0°	240	50	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 45°	240	50	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 90°	240	50	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 135°	240	50	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 180°	240	50	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 225°	240	50	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 270°	240	50	0	100	0	0.5 Cycle	(1)
0% UT for 0.5 Cycles @ 315°	240	50	0	100	0	0.5 Cycle	(1)
0% UT for 1 Cycle @ 0°	240	50	0	100	0	1 Cycle	(1)
70% UT for 25 Cycles @ 0°	240	50	70	30	168	25 Cycles	(1)
0% UT for 250 Cycles @ 0°	240	50	0	100	0	250 Cycles	(2)

Close Field Proximity Test Levels <i>Based on 0.3 m Separation Distance</i>			
Test Frequency (MHz)	Test Level (volts/meter)	Modulation (@ 50% duty cycle)	Communication Service (partial list)
385	27	18 Hz	TETRA 400
450	28	FM (5 kHz deviation)	GMRS/FRS
710, 745, 780	9	217 Hz	LTE
810, 870, 930	28	18 Hz	GSM 800
1720, 1845, 1970	28	217 Hz	GSM 1800
2450	28	217 Hz	RFID
5240, 5500, 5785	9	217 Hz	WLAN

GUIDANCE AND MANUFACTURER'S DECLARATION


The VuPad is intended for use in the electromagnetic environment specified below. The customer or user of the VuPad should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment Guidance
Radiated RF Emissions CISPR 11	Group 1	The VuPad uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Conducted RF Emissions CISPR 11	Class B	The VuPad is suitable for use in all establishments, including domestic, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonics IEC 61000-3-2	Class A	
Flicker IEC 61000-3-3	Complies	

The VuPad is intended for use in the electromagnetic environment specified below. The customer or user of the VuPad should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic Discharge [ESD] IEC 61000-4-2	± 8kV Contact ± 15kV Air	± 8kV Contact ± 15kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power Frequency [50/60 Hz] IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical commercial or hospital environment

The VuPad is intended for use in the electromagnetic environment specified below. The customer or user of the VuPad should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	10 V/m 150 kHz to 80 MHz	10 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the VuPad, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.2\sqrt{P}$ <p>$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 1.2\sqrt{P}$ 800 MHz to 2.7 GHz</p>
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m	<p>where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div>

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the VuPad is used exceeds the applicable RF compliance level above, the VuPad should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the VuPad.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10V/m.

WARNINGS AND CAUTIONS



WARNING: AN INSTRUCTION THAT DRAWS ATTENTION TO RISK OF INJURY OR DEATH

DANGER!

THE WATER PATH PROBE IS A TYPE B APPLIED PART AND PROVIDES A BASIC DEGREE OF PROTECTION AGAINST ELECTRIC SHOCK AND HAS DIRECT EARTH CONNECTION. WHEN USING THE WATER PATH PROBE IN UBM MODE (35 OR 50 MHZ), THE USER MUST TAKE EVERY PRECAUTION TO PREVENT THE TRANSDUCER FROM TOUCHING THE EYE. CONTACT BETWEEN A MOVING TRANSDUCER AND THE EYE CAN CAUSE SEVERE INJURY.

WARNINGS

Switching on a cold instrument near 0° Celsius may cause permanent damage. Allow the instrument to reach a normal room temperature for half a day in order to allow the internal elements to warm up and to avoid any thermal shock hazards when switched on. The cover will quickly reach room temperature, but not the internal circuitry.

To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth. Isolation from the supply mains may be achieved by disconnecting the main power cord from the supply outlet.

DO NOT disassemble, modify, or remodel the unit or accessories. This may cause unit damage, malfunction, electrical shock, fire, or personal injury.

DO NOT attempt to repair or service this instrument. Any repair or service to this instrument must be performed by experienced personnel who are trained by Sonomed Escalon. Attempts to repair or service the instrument may result in serious injury to the operator or patient.

Measurements should not be attempted when ocular integrity is questionable. The user needs to exhibit care in manipulating the measurement tip. Force should not be exerted against the eye.

Disconnect the AC POWER before cleaning the system.

If the device is used in the United States in the 240 volts mode, a Centered-Tapped Single-Phase 240V power supply must be used.

The transducers are fragile. Dropping or striking any probe can cause malfunctions; handle all probes with care. If a probe should be dropped, inspect it carefully for chips and cracks, and make a “test” scan on a known object. Damage to the front of the transducer will reduce efficiency, and may cause premature failure of the electronics or may cause damage to the cornea.

DO NOT USE PROBES IF TIP IS DAMAGED. ALWAYS EXAMINE PRIOR TO USE TO ENSURE PROBE INTEGRITY

This device is not intended for fetal use.

Never autoclave a transducer or expose it to high heat.

Do not attempt to connect the device to any accessories or supplemental equipment other than that provided by Sonomed Escalon. This could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation. Equipment connected to the system must be IEC 60601-1 or IEC 60950 compliant. Additionally, do not load any additional software onto the system without prior authorization from Sonomed Escalon (doing so may void warranty). Connection of system with a network is responsibility of user, including assurance of data integrity and related protections.

This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. Harmful interference to other devices can be determined by turning this equipment on or off. Try to correct the interference using one or more of the following: Reorient or relocate the receiving device, increase the separation between the equipment, connect the equipment to an outlet on a different circuit from that to which the other devices(s) are connected and or/ consult the factory field service technician for help.

Portable radio frequency communications equipment (including antenna cables and external cables) should be used no closer than 12 inches (30 cm) to any part of the VuPad, including cables supplied by Sonomed Escalon. Otherwise, degradation of the performance of this device could result.

Modifications to this instrument are not allowed. This may cause unit damage, malfunction, electrical shock, fire, or personal and/or patient injury.

Do not use the device together with HF surgical equipment. HF surgical equipment may be damaged, which may result in fire.

CAUTIONS

PATIENT DATA IS NOT SAVED IF UNIT IS TURNED OFF BEFORE SAVING SCANS

Data will be saved under the same patient name until another has been selected.

Care should be taken by user to minimize exposure of patient to ultrasound energy by keeping examinations as short as possible.

Position such that console is well ventilated with easy access to disconnect power cords, and do not block the power adapter in the event that disconnecting should be necessary.

Probes must be connected or disconnected only when the unit is switched OFF.

The unit should not be connected to a Multiple Power Switch Outlet (MPSO) which is also used to provide power to devices not intended to be used in the patient environment. Doing so may compromise electrical safety of the device.

Do not place the unit near heat sources such as a heater or operate in the presence of flammable anesthetics.

In order to prevent patient-to-patient transfer of infection, after each use disinfect the measurement tip following accepted clinical procedures. Refer to the Maintenance, Care and Service section regarding the use of disinfectants and for probe cleaning instructions.

Dispose of all products in accordance with local and national regulations and codes.

The VuPad conforms to the emissions and immunity requirements IEC 60601-1-2:2014, Conducted Emissions, Group 1, Class B.

Essential performance of the VuPad may be lost if the unit is adversely exposed to external electromagnetic disturbances resulting in loss of patient data. When tested for electromagnetic disturbances, the VuPad did not exhibit malfunction or degradation of performance when subjected to power frequency magnetic fields of 50 Hz and 60 Hz, but it is recommended that use in close proximity to or stacking of other electronic devices should be avoided because it could result in improper operation. If such use is necessary, the VuPad and other equipment should be observed to verify normal operations.

In the event adverse external electromagnetic disturbances causes the VuPad to lockup, the unit may require a system reboot by restarting.

Third Party Equipment: The use of third-party equipment, cables or accessories, not made or authorized by Sonomed Escalon, invalidates the warranty of the unit, and adversely affect the unit's safe operation.

CYBER SECURITY RECOMMENDATIONS

VuPad systems incorporate Microsoft 10 as the software Operating System, and, as such, makes available to Operators the full scope of Microsoft Windows 10 security features as defense against cyber security threats. Failure to maintain cyber security could result in compromised device functionality, loss of data availability or integrity, or exposure of other connected devices or networks to security threats.

Sonomed Escalon recommends the following minimum procedures be followed in order to maintain a basic level of cyber security:

- ① **Utilize Device Only for Intended Use.** Limit or prohibit use of device for any purpose other than ophthalmic ultrasound, including internet browsing and email, to limit potential exposure to cyber security risks.
- ② **Verify Windows 10 Firewall is Enabled.** Devices are shipped from the factory with the Windows Firewall on by default. To make sure it hasn't been turned off, follow these steps:
 1. Open Windows Firewall by clicking the **Start** button and then clicking the **Search** icon. In the search box, type **firewall**, and then click **Windows Firewall**.
 2. In left pane, click **Turn Windows Firewall On or Off**. If prompted for an administrator password or confirmation, type password or provide confirmation.
 3. Below each network location type, click **Turn On Windows Firewall**, and then click **OK**. It is recommended that the firewall be turned on for all network location types.

IMPORTANT: If device is connected to a network, ensure that the device is placed behind a strong network firewall.

- ③ **Verify Automatic Updating for Windows 10 Operating System Enabled.** With automatic updating, the Operator doesn't have to search for updates online or worry that critical fixes or device drivers for Windows might be missing from the system. Windows update automatically installs important updates as they become available. The Operator can also set Windows Update to install recommended updates automatically or to inform the Operator that they're available. The Operator can also choose whether to turn on Microsoft Update, which delivers updates for other Microsoft products. Optional updates, such as language packs and updates from Microsoft Update, aren't installed automatically. Windows Update won't add any apps to the system without prompting for permission.

To turn on automatic updating:

1. Open Windows Update by swiping in from the right edge of the screen (or, if using a mouse, pointing to the lower-right corner of the screen and moving the mouse pointer up), tapping or clicking **Settings**, tapping or clicking **Control Panel**, and then tapping or clicking **Windows Update**.
2. Tap or click **Change Settings**.
3. Under **Important Updates**, choose the option that you want.
4. Under **Recommended Updates**, select the **Give me recommended updates the same way I received important updates** check box, then click Apply.

IMPORTANT: In order for automatic Windows Update to function, the device must be continuously connected to the Internet. If the device is not connected, Windows updates will need to be performed manually. To do so, regularly go to Windows Update per step 1 above and then click **Check and Install Updates** button.

- ④ **Install Windows 10 Compatible Anti-Virus Program.** The Operator should utilize an antivirus and antimalware program, and keep it current by regularly downloading updates from the program manufacturer's website. Many of these programs update automatically and can help protect the system from spyware and malicious software.

- ⑤ **Enable Windows 10 BitLocker.** The Operator can use BitLocker Drive Encryption to help protect files on the entirety of the drive. BitLocker can help block hackers from accessing the system files they rely on to access sensitive data, or from accessing a disk drive by physically removing it from the system and installing it in a different one. New files are automatically encrypted when added to the disk drive that used BitLocker. However, if these files are copied to another drive or a different PC, they're automatically decrypted. BitLocker can encrypt the drive Windows is installed on (the operating system drive) as well as fixed data drives (such as internal hard drives). The Operator can also use BitLocker To Go to help protect all files stored on a removable data drive (such as an external hard drive or USB flash drive).

For more information on cyber security, please consult your IT support staff and/or visit Microsoft security www.microsoft.com/security.

WARRANTY

Sonomed Escalon warrants its products are free of defects of labor and material for two (2) years for consoles, one (1) year for probes and cables, and one (1) year for associated computer components such as monitors, keyboards, and mice.

The following items are not covered:

Physical damage to the console or probes due to misuse or shock.

Damage or data loss due to power failures or fluctuations. The use of a line-interactive UPS is recommended to avoid this type of failure.

Loss or corruption of data or software due to user error or the installation or use of any third-party hardware or software.

Damage to transducers caused by autoclaving or exposure to excessive heat.

Repairs not covered by warranty will be invoiced on the basis of parts and labor. At Sonomed Escalon's discretion, the damaged component may be exchanged at a flat rate.

Servicing of the unit may only be performed by Technicians certified by Sonomed Escalon. For additional information regarding system repair, maintenance, or exchange please contact US:

Sonomed Escalon
1979 Marcus Avenue, C105
Lake Success, NY 11042 USA
Tel: 800-227-1285
Fax: 516-354-5902
www.sonomedescalon.com

Any serious injury or incident occurring as a result of ophthalmic ultrasound use should be reported to Sonomed Escalon immediately. Patients and users should report serious incidents to the appropriate regulatory authorities.










NETWORKING

Sonomed Escalon does not provide support for the operation of this product in a network environment. Connection to and operation on any network is entirely the responsibility of the user. Where installation or use of any network hardware or software interferes with the normal operation of this Sonomed Escalon-supplied product, that product must be returned to normal operation at the user's expense. When the connection of this product to, or installation of Sonomed Escalon supplied software on, a network interferes with the operation of the network, the product must be removed from the network; alternatively, the problem may be resolved by the user in cooperation with the network owner, at their expense.

THIRD-PARTY SOFTWARE

Sonomed Escalon does not provide support for the use or installation of any software obtained from a third party on its products, including, but not limited to, operating system upgrades and device drivers. When software not supplied by Sonomed Escalon interferes with the operation of the system, the product will be returned to its original condition at the user's expense. Sonomed Escalon may occasionally furnish to users software not directly related to the functioning of its products. Such software is supplied as is, without warranty of any kind, and the availability of support for such software is at Sonomed Escalon's sole discretion.

SYMBOLS

	USB Port
	Do Not Dispose of Equipment in Normal Waste Stream
	CE Mark; Device complies with (EU) Medical Device Directive 93/42/EEC The number beside the symbol is the identification number of the Notified Body that certified the quality system under the Medical Device Directive
	Type BF Applied Part (B Scan Probe / A Scan Probe / Pachymeter Probe))
	Type B Applied Part (Water Path Probe)
	Warnings and Cautions, Read Accompanying Documents
IPX7	Rating applies to ultrasound probes only. An IPX7 rating means the probes are protected against effects of temporary immersion in water. Do not immerse probe electrical connectors.
IPX1	Rating applies to VuPad base unit only. An IPX1 rating means the device enclosure is resistant to the ingress of water at a flow rate of 1 mm/min for 10 minutes.
IPX6	Rating applies to foot pedal only. An IPX6 rating means the foot pedal can resist heavy sprays of high-pressure water.
	Class II Equipment per IEC EN 60601-1
	ETL Listing Mark
	Refer to Operator's Manual